

VTPC catalog



Communication software

Communication software

Software for communication creation

Purpose of use:

The software is suitable for individuals with speech, language, and communication impairments.

It provides the opportunity to convert text and/or symbols into spoken language, as well as allows the user to utilize the most appropriate access method for the software (touch, eye control, computer mouse/button).

Description:

- The software operates on a device with the Windows 10, Windows 11 operating system installed.
- Communication software can be used together with the Tobii Dynavox PCEye eye control device.
- The software includes various categories related to different vital needs of a person, daily activities, social needs, such as: emotions, food, clothing, personal needs, activities, and many other categories.
- Individual settings are available, tailored to the user's abilities.
- Environmental control using infrared radiation (TV, toys, adjustable beds, lighting, air conditioners, etc.).
- It is possible to make audio recordings using one's voice, as well as add text and symbols (images, photographs, etc.), according to the user's abilities and needs.
- Personalization options for text and symbol selections.
- The software includes synthesized speech, with the option to change the voice type (female, male, adult, child).

Additional equipment:

Software add-ons are available as needed - Latvian language license. (a physical person who uses the Latvian language in daily communication is provided the opportunity to set the voice function in the software to Latvian).

Important:

A personal computer is required for obtaining and using the software.

The device's software is compatible with the Windows operating system.

- ID number : 222112-0
- ISO code : 22 21 12
- Manufacturer : EU



Speech board with options

Speech board with options

Device for speech therapy

Purpose of use:


The alternative communication technical aid supports individuals with speech, language, communication, and movement disorders.

The device allows users to record short audio messages using their voice, which are enhanced with visual support - an image, pictogram, or photograph.

Individuals can play back audio recordings in various everyday situations (*at home, at school, in kindergarten, during sessions with specialists, etc.*), to express their wishes, needs, engage in dialogue, provide answers and comments, learn new skills, and facilitate communication with others.

Technical aid description:

- Operates at five easily switchable levels;
 - Four choice boxes and two basic recording boxes;
 - Total number of recordings - 20;
 - Total recording duration - 4.6 min;
 - Built-in handle for easy portability;
 - Volume control;
 - Built-in keypad for easier use;
 - Quick audio recording and deletion;
 - Device operates on two AA batteries.
 - Includes a CD for creating printouts.
-
- ID number : 222109-2
 - ISO code : 22 21 09
 - Manufacturer : USA



Speech board with options

Speech board with options

Device for speech therapy

- The technical aid is designed for individuals with speech and language disorders, as well as for individuals with movement and coordination difficulties, to engage in communication.
- The device allows recording short informative audio materials, which are enhanced with visual support in the form of images, pictograms, or photographs, and also allows patients to integrate into society,

daily activities, academic environments, and professional activities.

1. Operates at five easily switchable levels.
2. Nine choice boxes and three basic recording boxes.
3. The total number of messages is no less than 45.
4. The duration of one message is 8 seconds at any level.
5. The total recording duration is no less than 8.5 minutes.
6. Device with integrated volume control.
7. Handle for easy transportation.
8. Additional CD for creating prints that include 1000 images (pictograms), as well as allows adding your own images or photographs.
9. With the option to add text (in any language), change font, color, sizes, and also save in memory.

- ID number : 222109-2
- ISO code : 222109
- Manufacturer : EU



Speech board with options

Speech board with options

Device for speech therapy

Purpose of use:

The alternative communication technical aid supports individuals with speech, language, communication, and mobility impairments.

The device allows users to record short audio messages using their voice, which are enhanced with visual support - an image, pictogram, or photograph.

Individuals have the opportunity to play audio recordings in various everyday situations (*at home, school, kindergarten, sessions with specialists, etc.*), to express their wishes, needs, engage in dialogue, provide answers and comments, learn new skills, and facilitate communication with others.

Technical aid description:

- Operates at five easily switchable levels;
 - 20 choice boxes and 5 basic recording boxes;
 - Total number of recordings - 100;
 - Total recording duration 15 minutes;
 - Built-in handle for easy portability;
 - Volume control;
 - Built-in keypad for easier use;
 - Quick audio recording and deletion;
 - The device operates on two AA batteries.
 - Includes a CD for creating prints.
-
- ID number : 222109-2
 - ISO code : 22 21 09
 - Manufacturer : USA



Alternative communication device with eye tracking (Tobii Dynavox PCEye)

Alternative communication device with eye tracking (Tobii Dynavox PCEye)

Creating communication using eye gaze.

Purpose of use:

An alternative communication device that allows individuals with severe mobility impairments to control a computer and access all its functions using eye gaze.

It replaces the conventional keyboard and mouse, provides access to the internet, social networks, allows writing emails and creating documents, and enables adjusting the pace according to the user's abilities.

Technical aid description:

- The device is easily connectable to a computer that has *Windows 10, Windows 11 operating system installed.*
- The device can also be used by individuals with corrective glasses and contact lenses, both binocular and monocular.
- This device is suitable for individuals with severe mobility restrictions, for example, in cases of conditions such as – cerebral palsy, ALS, Rett syndrome, various degenerative progressive diseases, etc., resulting in severe motor function limitations or loss.
- For individuals who need symbol and/or text support for communication learning and creation, the device can be used together with communication software such as: *Communicator 5, TD Snap* etc.

Important:

A personal computer is required to receive and use the device.

The device's software is compatible with the Windows operating system.

- ID number : 223621-0
- ISO code : 22 36 21
- Manufacturer : EU



Communicator with Illuminated Choice Windows (Go Talk Express 32)

Communicator with Illuminated Choice Windows (Go Talk Express 32)

Device for speech therapy

- The technical aid is designed for individuals with speech and language impairments to engage in communication.
- The device allows recording short informative audio materials, which are enhanced with visual support in the form of images, and also allows patients to integrate into society, daily activities, academic environments, and professional activities.
 1. The device operates in two modes:
 - standard – plays one message,
 - express – plays multiple messages simultaneously, allowing for the playback of phrases, sentences, and narratives.
 2. For visual reinforcement – 4 LED lights around each of the message boxes.
 3. The device operates at five levels, has integrated volume control, and if necessary, levels can be locked.
 4. Handle for easy transportation.
 5. Includes a CD for creating printouts, which contains 4500 images (pictograms), and allows adding additional images or photographs.
- ID number : 222109-5
- ISO code : 222109
- Manufacturer : EU



High-level communicator

High-level communicator

- The device is intended for individuals with speech and language impairments.
- The device's message boxes are small (1.5x1.5cm), designed for individuals with well-developed fine motor skills.
- The device provides the opportunity to create semantically quality speech for individuals with speech and language impairments. Thanks to the device's fast operation, it is possible to form sentences, small sequences of sentences without long pauses between words.
- The multi-level communication device operates by playing pre-recorded messages, supplemented with images/pictograms for visual representation.
- ID number : 222109-3
- ISO code : 222109
- Manufacturer : EU



Communication Software

Communication Software

Software for Communication

Purpose:

Alternative communication software suitable for individuals with speech, language, or communication impairments.

The program allows users to express wishes, needs, and communicate with others. Selected symbols are played from the device, generating speech.

The program is suitable for individuals with autism spectrum disorders, cerebral palsy, Down syndrome, intellectual disabilities, aphasia, and other impairments.

Description:

- Adjustment of the number of message boxes;
- Creation of individual spreadsheets according to the user's needs and abilities;
- Addition of personal images;
- Input options using touch, buttons, mouse, scanning, or eye gaze control;
- The software operates on a device with the Windows operating system installed.

Additional Equipment:

As needed, additional software equipment is available - Latvian language license. *(A physical person who uses the Latvian language in daily communication is provided with the option to set the voice function of the software in Latvian).*

Important:

A personal computer is required for receiving and using the software.

The device's software is compatible with the Windows operating system.

- ID number : 222112-0
- ISO code : 22 21 12
- Manufacturer : EU



Programmable communicator (Tech/Scan 8 Plus)

Programmable communicator (Tech/Scan 8 Plus)

The assistive device is designed for individuals with speech and language impairments.

- The device's configuration allows not only for the classic recording and playback of audio material but also for modulating the playback speed.
- The multi-level communication device operates by playing pre-recorded messages with visualization options.
- 8 message playback slots, 8 levels, total number of messages 64.
- ID number : 222109-4
- ISO code : 222109
- Manufacturer : EU



Multilevel Speech Device

Multilevel Speech Device

Device for speech therapy - for making two choices.

The device is intended for individuals with speech, language, coordination, and movement disorders to engage in communication.

- The device has two easily pressable, noticeable, colorful buttons, which the patient can use to express their choice.
- The device operates on three levels, with a total of no less than six messages.
- The total recording time is no less than 6 minutes, with at least 2 minutes for each message.
- The device has integrated image/pictogram storage. (see att.)

- ID number : 222109-0
- ISO code : 222109
- Manufacturer : EU



Communicator for playing recorded messages

Communicator for playing recorded messages

The device for speech therapy is intended for individuals with speech and language impairments.

- The device comes with four interchangeable message boxes of varying button counts.
- The equipment configuration allows for individual customization of communication fields according to each patient's abilities and needs, as well as changing this choice during use as necessary.
- Number of message boxes - 1, 2, 4, and 8.
- The device operates at 4 levels, corresponding to the selected number of messages to be played.
- ID number : 222109-1
- ISO code : 222109
- Manufacturer : EU